

There are three models available :

I. Bioptron 1 -With circular therapy area of 5 cm for use for smaller problem area, more for a home kit :

II. Bioptron 2 - With circular area of 15 cm commonly used in the clinics and hospitals. It comes along with a trolley or a wall mount support.

III. Bioptron 3 - Equipment for the whole body therapy with therapeutic area covering 30 x 120 cm, specially for use in the major hospitals, in the Pain Departments, Burn Departments, in conditions like psoriasis, neurodermatitis and SLE conditions.

OXY-spray is applied in combination with Bioptron Light Therapy. This spray contains pure water enriched with medical oxygen. It provides the necessary oxygen content to the area of therapy. When used with the Bioptron Light Therapy, it gets absorbed readily and act accordingly in building immunity and repairing of the skin on the treated area

**BIOPTRON**

*Certified*

**Medical Device**

*by*

THE INTERNATIONAL CERTIFICATION  
NETWORK

ISO 9001:1994/EN 46001:1996  
Reg. No. CH-14356-02

**IQNet and SQS**

*IQNet is represented in the USA by AFAQ AIB Vingotte  
International CISQ, DQS, KEMA, NSAI and QMI*

BIOPTRON AG  
Sihleggstrasse 23  
CH-8832 Willerau  
☎ +41 43 868 28 00  
☎ +41 43 868 28 99  
✉ [light@bioptron.com](mailto:light@bioptron.com)



**EC REP**

Zepter Handelsgesellschaft mbH  
Pummenstrasse 34  
A-4020 Linz

**zepter**  
INTERNATIONAL  
zepter

**BIOPTRON**  
[www.bioptron.com](http://www.bioptron.com)

Local contact

*Bains Physio*

**zepter**  
**BIOPTRON**   
LIGHT THERAPY SYSTEMS



ENERGY WITH  
LIGHT THERAPY